

medatro[®] medical trolley

B28

medical industry mobility solution

Safety • Durability • Mobility • Customized LOGO

Dimension: φ600*890mm

Base size: φ600*70mm

Column size: 78*100*810mm

Platform: 275*200*5mm

Infusion pump stand: 275*200*5mm

Basket: 330*210*135mm

Handle: 230*175mm

Rob of handle: φ30*222mm

Cable hanger: 204*65*34mm

Castor: 5pcs*4inch (with brake)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Net weight: 19.4kg

Carton dimension: 900*570*210mm

Gross weight: 22kg

