

medatro® medical trolley**K10**

medical industry mobility solution

Safety · Durability · Mobility · Customized LOGO

Dimension: 600*550*850mm

Base size: 600*550*165mm

Column size: 70*135*850mm

Platform: 275*200*5mm

Infusion pump bracket: $\varnothing 25 \times 700$ mm

Handle: $\varnothing 22 \times 450 \times 160$ mm (red)

Installing plate for handle: 130*100*8mm

Castor: 4pcs*4inch (with brake)

Load Capacity: 30kg

Max push speed: 2m/s

Max tilt angle: 15°

Carton dimension: 805*635*300mm

Net weight: 20kg

Gross weight: 22.3kg

