

medatro<sup>®</sup> medical trolley

## G01

medical industry mobility solution

**Safety • Durability • Mobility • Customized LOGO**

**Dimension:** 510\*550\*1100mm

**Base size:** 510\*550\*30mm

**Column size:**  $\phi 60$ \*893\*2.5mm

**Platform:** 295\*240\*10mm

**Humidifier hanger:** 55\*40\*16mm

**Handle:** 170\*35\*45mm

**Rob of handle A:** 100\*50\*25mm

**Rob of handle B:** 100\*50\*25mm

**Fixed bar:** 440\*25\*10mm

**Castor:** 4pcs\*4inch (with brake)

**Load Capacity:** 30kg

**Max push speed:** 2m/s

**Max tilt angle:** 15°

**Net weight:** 15.4kg

**Carton dimension:** 900\*570\*210mm

**Gross weight:** 18kg

